



CHROMOPTIC

La Maison des Services - Avenue des 2 lacs
ZAC de Villejust - BP 40139
91942 VILLEBON COURTABOEUF Cedex



www.chromoptic.com



contact@chromoptic.com



01.69.31.41.65

SOLUTION D'ÉTALONNAGE

6610030100 (AGILENT TECHNOLOGIES)



Solution d'étalonnage des longueurs d'ondes des ICP-OES

Contient 5mg/ml de Al, As, Ba, Cd, Co, Cr, Cu, Mn, Mo, Ni, Pb, Se, Sr et Zn

Contient 50mg/ml de K

(Aluminium nitrate nonahydrate ICP grade CAS# 7784-27-2; Barium nitrate ICP grade CAS# 10022-31-8; Potassium Nitrate CAS#10042-76-9; Strontium Nitrate CAS# 10042-76-9)

Solvant: Acide Nitrique HNO₃ à 5%

Conditionnée en bouteille de 500ml