



## CHROMOPTIC

La Maison des Services - Avenue des 2 lacs  
ZAC de Villejust - BP 40139  
91942 VILLEBON COURTABOEUF Cedex



[www.chromoptic.com](http://www.chromoptic.com)



[contact@chromoptic.com](mailto:contact@chromoptic.com)



01.69.31.41.65

### COLONNES DE GARDE POROSHELL

821725-911 (AGILENT TECHNOLOGIES)



Colonne InfinityLab Poroshell 120

Phase inverse EC-C18

Longueur: 5mm

DI: 2.1mm

Taille des particules: 2.7µm, partiellement poreuses

Porosité: 120Å

Pression max: 600 Bar

Taux de carbone: 10%

Postsilanisée

Température max: 40°C à pH 6-8 et 60°C à pH 2-6

Limites de poids moléculaire: 0 à 3000 Da

Solvant de stockage: acétonitrile / eau

Gamme de pH: 2-8

Conditionnées par 3