






## CHROMOPTIC

La Maison des Services - Avenue des 2 lacs  
ZAC de Villejust - BP 40139  
91942 VILLEBON COURTABOEUF Cedex

 [www.chromoptic.com](http://www.chromoptic.com)  
 [contact@chromoptic.com](mailto:contact@chromoptic.com)  
 01.69.31.41.65

### COLONNES METAGARD PURSUIT PFP

A3051MG2 (AGILENT TECHNOLOGIES)



Colonne HPLC Pursuit

Phase : PFP

DI: 2mm

Granulométrie: 3 $\mu$ m, partiellement poreuses

Volume de pores: 11ml/g

Couverture ligand: 3.4 $\mu$ mol/m<sup>2</sup>

Porosité: 200Å

Surface: 200m<sup>2</sup>/g

Taux de carbone: 6.3%

Postsilanisée

Gamme de pH: 1.5 -10

Conditionnées par 3