

Separations | HT8 GC column

More than five decades of innovative phase technologies and unique fused silica production capabilities, together provide end-to-end separation solutions for all applications.

Features and benefits

- Preferred column for polychlorinated biphenyl (PCB) compounds.
- Separates PCB's on ortho ring substitution as well as boiling point.



Recommended applications

- PCB congener analyses.
- Nitro-substituted aromatics.
- Polynuclear aromatic hydrocarbons.
- Pesticides/herbicides.

Product specifications

8% Phenyl Polycarborane-siloxane. Suitable replacement for: no equivalents, unique ultra high temperature column.

Film thickness	Operating temperature
0.25 µm	-20°C to 320/350°C



Environment

For more information about this product visit www.trajanscimed.com or contact techsupport@trajanscimed.com

Contactez-nous pour plus d'informations:

Région Sud-Est: **Hervé BRETON**

h.breton@chromoptic.com 06.75.55.98.45

Région Ile-de-France Centre: **Céline LAURENCE**

claurence@chromoptic.com 07.85.27.58.57

Région Nord et Est: **Valérie FILLEBEEN**

vfillebeen@chromoptic.com 06.76.68.29.65

Région Grand Ouest: **Marjolaine LEBON**

mlebon@chromoptic.com 06.36.77.11.25