



YMC-Triart C18




versatile applications
first choice for
method development
pH-stable 1 - 12
also 100% aqueous
eluents

YMC-Triart C18 ExRS



pH-stable 1 - 12
hydrophobic substances
positional isomers


YMC-Triart C8




alternative to C18
short retention time
pH-stable 1 - 12

Easy YMC-Triart Column Selection

for Method Development




YMC-Triart Phenyl



aromatic compounds
(π -electron acceptor)
conjugated systems
also 100% aqueous
eluents

YMC-Triart PFP



aromatic compounds
(π -electron donor)
cis-trans isomers
polar halogenated
compounds
also 100% aqueous
eluents

YMC
EUROPE GMBH

Schöttmannshof 19
D-46539 Dinslaken
Phone: +49 (0)2064 427-0
Fax: +49 (0)2064 427-222
E-Mail: info@ymc.de


www.ymc.de




CHROMOPTIC

CHROMOPTIC SAS · La Maison des Services
ZAC de VILLEJUST · Avenue des 2 Lacs
91971 COURTABOEUF CEDEX
Phone: +33 69 31 41 65
E-Mail: contact@chromoptic.com


www.chromoptic.com




YMC-Pack Pro C18



YMC-Pack Pro C18 RS




Hydrosphere C18




Easy ProFamily Column Selection


for Method Development



YMC-Pack Pro C8



YMC-Pack Pro C4



YMC
EUROPE GMBH

Schöttmannshof 19
D-46539 Dinslaken
Phone: +49 (0)2064 427-0
Fax: +49 (0)2064 427-222
E-Mail: info@ymc.de

www.ymc.de

CHROMOPTIC

CHROMOPTIC SAS · La Maison des Services
ZAC de VILLEJUST · Avenue des 2 Lacs
91971 COURTABOEUF CEDEX
Phone: +33 69 31 41 65
E-Mail: contact@chromoptic.com

www.chromoptic.com

